

Act ecologically:
The TOA Mobile Toilet
The Biological Dry Toilet for every need



- requires minimal space
- offers same comfort as usual flush toilets
- needs no chemicals or water
- creates hygienic sanitation everywhere, where no toilet is
- is easy to handle and clean
- offers good conditions for composting
- bucket manufactured of frost-proof polyethylene-plastic
- solid design



Modell TOA Standard



Modell TOA Mini



Modell TOA Maxi

The TOA Mobil Toilet

There are three models of the TOA Mobile Toilet: The **TOA Standard Mobile Toilet** consists of a bucket, a tight cover and a toilet seat. The **TOA Maxi, TOA Luxus** as well as the **TOA Mini Mobile Toilet** is additionally equipped with a housing, containing the bucket. The toilet seats can be opened for example to take out the bucket. The toilet seat of the Standard model is fastened on top of the bucket and can be easily removed for emptying.

1. Function

Before the first use, a bedding of absorbing and structural material such as bark chips, wood clippings, earth or shredded material or a mixture of all is put into the bucket. Take as much bedding as fluid is expected until the next emptying, so the solid waste is not in contact with the fluid level. After each defecation the toilet paper should be covered with a handful of bedding material. This way, waste is hidden from view and odours are largely bound over a period of at least 1-2 days. For transport, the bucket can be sealed up with a tight cover.

2. Range of applications

The TOA Mobile Toilet is designed for use in all purposes, cases of necessity as well as areas of spare time activities and every other purpose, where easy and reasonable solutions are demanded and bedding material is provided. Due to easy transport, composting is independent from toilet location.

3. Building conditions

The TOA Mobile Toilet needs minimal space for use. Toilet rooms should be aerated as any other toilet room.

4. Plumbing

The TOA Mobile Toilet is independent from water and electricity.

5. Maintenance

To keep the bucket clean very easily, we offer bags for both models, which can be emptied onto the compost as well. The bags decompose within six months. Larger capacities can be served with another bucket, which can be exchanged when the first bucket is filled.

<u>Dimensions:</u>	Standard	42 cm wide, 43 cm deep, 42 cm high, volume 30 litres
	Mini	35 cm wide, 44 cm deep, 43 cm high, volume 12 litres
	Luxus	39 cm wide, 50 cm deep, 46 cm high, volume 25 litres
	Maxi	45 cm wide, 54 cm deep, 52 cm high, volume 35 litres

<u>Material:</u>	Standard	bucket + cover of polyethylene-plastic, beige seat of pvc- plastic, white
	Mini	housing (grey) + bucket (black) of polyethylene-plastic, seat of pvc-plastic, ruby coloured
	Luxus	bucket + cover of polyethylene-plastic, white seat + housing of abs-plastic, white
	Maxi	bucket + housing of polyethylene-plastic, white seat + cover of polypropylene- plastic

BERGER BIOTECHNIK® GmbH

Bogenstr. 17
D - 20144 Hamburg
Tel.: +49 (0) 40 / 439 78 75
Fax: +49 (0) 40 / 43 78 48
eMail: info@berger-biotechnik.de
Internet: www.berger-biotechnik.com